

---

# An introduction to Particle Image Velocimetry (PIV)

## *Principles and Applications*

---

***Fabio Cozzi***

Laboratorio di Combustione



Dip. di Energia, Politecnico di Milano



*April 2014*

# Imaging Systems


Imaging systems are non intrusive optical techniques for research and diagnostics into flow, turbulence, micro fluidics, spray atomisation and combustion processes.



# Quantitative Flow Visualization



- Quantitative flow visualization requires the fluid must be seeded with small particle tracers able to follow the instantaneous changes of the flow velocity.
- In order to obtain robust unbiased measurements over the flow domain, it is important that tracer particles are homogeneously distributed within the observed flow region.



Visualization vs. Measurement

# What is PIV?



- PIV stands for *Particle Image Velocimetry*
- It is an optical diagnostic technique used to measure the velocity field in gas/liquid flows
- Actually the displacement of particles scattered in the flow is measured!

# The Principle of PIV

The principle of PIV is based on the measurement of the displacement of small tracer particles that are carried by the fluid during a short time interval.

The tracer particles are sufficiently small that they accurately follow the fluid motion and do not alter the fluid properties or flow characteristics

PIV can be applied to virtually any kind of flow, as long as the fluid is *transparent* and a suitable optical access to the flow is available





Fig. 1. Schematic of the PIV method.

# The Principle of PIV




The basic principle is to take two snapshots of the flowfield at a known time interval, and evaluating the velocity

Once a sequence of two light pulses is recorded, the images are divided into small subsections called interrogation areas ( $D_I$ ).



The velocity vectors are statistically derived from each interrogation area by measuring the average displacement of particles between two light pulses:



$$\vec{V} = \frac{\Delta \vec{s}}{\Delta t}$$

# Why use PIV?



## ADVANTAGES

- ↪ non-intrusive (besides particles)
- ↪ 2 velocity components are readily measured (in a plane!)
- ↪ rapid mean flow measurement
- ↪ allow flow structures to be visualized
- ↪ spatial correlations
- ↪ extension to
  - ✓ three velocity components
  - ✓ volume mapping


## LIMITATIONS

- ↪ needs optical access
- ↪ max area of interest about 0.5 m x 0.5 m
- ↪ reflections at walls
- ↪ safety
- ↪ relatively expensive



# Components of PIV System

- Light source
  - Laser ....
  - Light sheet formations optics
- Image recording system
  - Camera (film, *CCD*, *CMOS*)
  - Optics (objective + filter)
  - Intensifier (not always needed)
  - Frame grabber (not always needed)
- Synchronization Unit
- PIV Image Recording and Processing System
  - PC
  - Storage systems
  - Software
- Seeding system





# SEEDING (introduction)




- ❑ A very important issue for obtaining accurate PIV measurements is appropriate seeding of the flow with tracer particles.
- ❑ To closely follow the flow the particles should be as small as possible, but on the other hand they may not be too small, because then they will not scatter enough light, and hence produce too weak images.
- ❑ For liquid flows a particle diameter of  $d_p = 10\text{-}20\ \mu\text{m}$  is a usual compromise.
- ❑ In gas flows the particles have to be smaller to follow the flow, because of the lower density of gas compared to liquid. In gas flows a usual compromise of particle size is  $d_p = 1\text{-}5\ \mu\text{m}$ .
- ❑ **It is desirable that seeding particles be**  
non-toxic, non-corrosive, non-abrasive, non-volatile and chemically inert.

# SEEDING (particle motion in flows)



- At low Reynolds number the motion of a spherical particle is describe by BBO equation



$a$  = particles radius

$v_p$  = particles velocities

$U$  = gas velocities

$V = v_p - U$

$$\begin{aligned}
 & \frac{4}{3}\pi a^3 \rho_p \frac{dv}{dt} \\
 = & \underbrace{\frac{4}{3}\pi a^3 \rho \frac{DV}{Dt}}_{\text{non-inertial force}} + \underbrace{\frac{4}{3}\pi a^3 (\rho_p - \rho_f) g}_{\text{net body force}} \\
 & - \underbrace{6\pi\mu a \left[ (v_p - U) - \frac{1}{6}a^2 \nabla^2 U \right]}_{\text{quasi-steady viscous force}} \\
 & - \underbrace{6\pi\mu a^2 \int_0^t \frac{d}{d\tau} \left[ (v_p - U) - \frac{1}{6}a^2 \nabla^2 U \right] \frac{d\tau}{\sqrt{\nu(t-\tau)}}}_{\text{time history (Basset) force}} \\
 & - \underbrace{\frac{2}{3}\pi a^3 \frac{d}{dt} \left[ (v_p - U) - \frac{1}{10}a^2 \nabla^2 U \right]}_{\text{added mass force}} + \underbrace{\frac{1}{10}\pi a^3 \nabla^2 U}_{\text{lift force}},
 \end{aligned}$$

Basset–Boussinesq–Oseen (BBO) equation

# SEEDING (particle motion in flows)



According to previous approximations a small particle moving in a fluid is described by the equation of motion

$$\frac{4}{3}\pi a^3 \rho_p \frac{dv_p}{dt} = \frac{C_D}{2} \pi a^2 \rho_f V^2 \quad \text{Re}_p = \frac{2a\rho_f V}{\mu}$$

$\mu$ : fluid viscosity,  $C_D$ : drag coefficient,  $V$  relative gas velocity

In the Stokes regime ( $\text{Re}_p < 1$ )  $C_D = 24/\text{Re}_p$  and if  $\rho_p \gg \rho_f$

The particle has a characteristic response time:  $\tau = \frac{1}{18} \frac{\rho_p d_p^2}{\mu}$

and assumes the local fluid velocity exponentially.

In swirling flows the seeding particles are forced outwards in a manner that depends critically on the particle diameter, particle mass and vortex circulation. In this case the particle characteristic time should be less than the vortex time scale.

# SEEDING (particle motion in flows)



## ■ Example

100  $\mu\text{m}$  oil particle in air  $(\rho_p=900 \text{ kg/m}^3)$   $\tau \sim 0.03 \text{ s}$

10  $\mu\text{m}$  oil particle in air  $(\rho_p=900 \text{ kg/m}^3)$   $\tau \sim 0.0003 \text{ s}$

10  $\mu\text{m}$   $\text{Al}_2\text{O}_3$  particle in water  $(\rho_p=3960 \text{ kg/m}^3)$   $\tau \sim 0.00003 \text{ s}$

# SEEDING (particle motion in flows)



$$\rho_p \gg \rho_f$$

$$\frac{U_p}{U_f} = \frac{1}{\sqrt{1 + \omega^2 \tau^2}}$$


$$\tau = \frac{1}{18} \frac{\rho_p d_p^2}{\mu}$$

**Table 2.** Particle response in turbulent flow ( $U_p/U_f = 0.99$ )

Particle	$\rho_p$ (kg m <sup>-3</sup> )	Gas (10 <sup>5</sup> Pa)	Density ratio $s$	Viscosity $\nu$ (m <sup>2</sup> s <sup>-1</sup> )	$f_c$ (kHz)	$Sk_c$	$d_p$ ( $\mu$ m)
TiO <sub>2</sub>	3500	Air (300 K)	2950	$1.50 \times 10^{-5}$	1	0.0295	1.44
					10		0.45
Al <sub>2</sub> O <sub>3</sub>	3970	Flame (1800 K)	20250	$3.00 \times 10^{-4}$	1	0.0113	2.46
					10		0.78
Glass	2600	Air (300 K)	2190	$1.50 \times 10^{-5}$	1	0.0342	1.67
					10		0.53
Olive oil	970	Air (220 K)	617	$1.45 \times 10^{-5}$	1	0.0645	3.09
					10		0.98
Microballoon	100	Air (300 K)	84.5	$1.50 \times 10^{-5}$	1	0.1742	8.50
					10		2.69

$$Sk = \left(\frac{\omega}{\nu}\right)^{1/2} d_p$$

$$\nu = \mu/\rho_f$$



# SEEDING



For PIV-LDV applications in isothermal flows the best solution is to filter the air and introduce artificially generated aerosol particles.

The injection of seeding particles can perturb the reacting flow through:

- Catalytic effects
- Heat exchange (cooling of the gas)
- Absorption of light (obscuration)
- Particles/vortices interactions (typically particles are smaller than the smallest turbulence scales)

These effects depend also on the particle number concentration.

Typical probe volume dimension:  $10^{-11} \text{ m}^3$

Max.  $N_p \approx 10^{11} \text{ m}^{-3}$  (single scattering, LDV)

Particle concentration determines the average Data Rate and frequency resolution of the turbulent velocity spectra.

# SEEDING



## ■ GAS FLOW

Material	particle diameter $d_p$ ( $\mu\text{m}$ )	Laser	Pulse energy, pulse time	Light sheet	
				$w$ (mm)	$t$ (mm)
TiO <sub>2</sub> ( $m = 2.6$ , $\rho = 3500 \text{ kg m}^{-3}$ )	<1	Nd:YAG	10 mJ, 20 ns	15	0.3
TiO <sub>2</sub> , ZrO <sub>2</sub>	0.7–1	Nd:YAG	110 mJ, 12 ns		
Al <sub>2</sub> O <sub>3</sub> ( $m = 1.76$ , $\rho = 3970 \text{ kg m}^{-3}$ )	0.3	Nd:YAG	400 mJ		0.2
	3	Nd:YAG	9 mJ, 6 ns		
	0.8	Ruby	20 ns	150	$\simeq 1$
Polycrystalline	30	Nd:YAG	135 mJ, 6 ns		
Glass	30	Ruby	30 mJ, 30 ns		
Oil smoke	1	Ruby	5 J		
Corn oil	1–2	Nd:YAG	100 mJ		
Oil	1–2	Nd:YAG	120 mJ		0.4
Olive oil ( $m = 1.47$ , $\rho = 970 \text{ kg m}^{-3}$ )	1.06	Nd:YAG	70 mJ, 16 ns	200	0.5

Fluid	Material	Diameter ( $\mu\text{m}$ )	Density ( $\text{kg/m}^3$ )
Air	DEHS	1–3	$10^3$
–	Glycol–water solution	1–3	$10^3$
–	Vegetable oil	1–3	$10^3$
–	TiO <sub>2</sub>	0.2–0.5	$1\text{--}4 \times 10^3$

Material	Density ( $\text{g/cm}^{-3}$ )	Index of refraction	Melting point ( $^{\circ}\text{C}$ )
NaCl	2.16	1.54	801
Al <sub>2</sub> O <sub>3</sub>	3.96	1.79	2015
TiO <sub>2</sub>	4.26	2.6–2.9	1750
SiC	3.2	2.6	2700
ZrO <sub>2</sub>	5.6	2.2	2980



# SEEDING





## ■ LIQUID FLOW

- ❑ The higher density/viscosity of liquids as compared to gases allows to use bigger particles
- ❑ Bigger particles require less laser power

Material	particle diameter $d_p$ ( $\mu\text{m}$ )	Laser	CW power or energy, time	Light sheet	
				$w$ (mm)	$t$ (mm)
TiO <sub>2</sub>	3	Nd:YAG			
Al <sub>2</sub> O <sub>3</sub>	9.5	Ruby	2 J, 30 ns	100	0.8
Conifer pollen ( $\rho = 1000 \text{ kg m}^{-3}$ )	50–60	Ar ion	1–2 W		
Polymer ( $\rho = 1030 \text{ kg m}^{-3}$ )	30	Ar ion	0.5–5 W		0.5
Phosphorescent polymer	80	Ar ion	5 W		1
Fluorescent	50	Nd:YAG			
	20	Cu vapour	45 W		1
Polystyrene ( $\rho = 1050 \text{ kg m}^{-3}$ )	500				
	15	Ruby	25 mJ, 20 ns		
Thermoplastic ( $\rho = 1020 \text{ kg m}^{-3}$ )	6	Nd:YAG		50	2
Reflective ( $\rho = 1010 \text{ kg m}^{-3}$ )	60	Ar ion	18 W		
	30	Ar ion	12–18 W	200	
Metallic coated	4	Ar ion	2 W		2
	14	Ar ion			1
Microspheres ( $\rho = 700 \text{ kg m}^{-3}$ )	<30	Ar ion			
H <sub>2</sub> bubbles		Ar ion	1 W		0.3

# SEEDING


## ■ Aerosol generation



**Fig. 5.65** Schematic of a Laskin nozzle unit for droplet seeding

# SEEDING

## ■ Aerosol generation




**Fig. 5.67a,b** Aerosol generation from powder. **(a)** Operation principle of a rotary brush seeder. **(b)** A cyclone aerosol generator

# SEEDING (light scattering)




- Laser light is scattered everywhere by the particle
  - The total amount of scattered light,  $I_0$ , is a function of:  
Refractive index, light polarization plane, particle diameter, particle shape .....



# SEEDING (light scattering)



- Computed light scattering cross section as a function of non-dimensional particle diameter




# SEEDING (light scattering)




## ■ LIGHT SCATTERING FROM SMALL PARTICLES


- Particles should scatter enough light to be recorded and to produce high contrast images
- Light scattered by small particles is a function of
  - Ratio of the refractive index of the particles to that of the surrounding medium scattering efficiency  $\uparrow$  as this ratio  $\uparrow$
  - Particle size, shape and orientation
  - Light polarization and observation angles (for PIV usually  $90^\circ$ )
- For  $d_p > \lambda$  Mie's scattering theory apply
  - Light scatter is not uniform! If  $q = \pi d_p / \lambda > 1$ , approximately  $q$  local maxima appear in the angular distribution over the range from  $0^\circ$  to  $180^\circ$
  - By increasing the particle size (i.e.  $\pi d_p / \lambda$ ) the ratio of forward to backward scatter intensity will increase rapidly.



$$d_p \approx 0.2\lambda$$



$$d_p \approx 1.0\lambda$$



$$d_p \approx 10\lambda$$

# Imaging of small particles

The image of a point source (scattering particle) does not appear as a point in the image plane of a camera, but forms a Fraunhofer diffraction pattern (Airy disk), even if it is imaged by a perfect aberration-free lens.





Figure 5. Diffraction effects during particle image recording.



Airy pattern of a small and a larger aperture diameter



# Imaging of small particles

The diameter,  $d_{diff}$ , of the Airy disk represents the smallest particle image that can be obtained by a given imaging configuration:

$$d_{diff} = 2.44 f_N (M + 1) \lambda \quad f_N = f / D$$


$M = z_o / Z_0$  magnification factor;

$f$  = focal length of the lens;  $D$  = aperture diameter

The image of a particle is given by the convolution of the Airy disk and the geometric image of the particle ( $d_g = M d_p$ ):

$$d_i = \sqrt{(M d_p)^2 + d_{diff}^2}$$

Normalized intensity distribution of the Airy pattern



# Imaging of small particles



Moreover, the particle-image diameter follows the previous equation for particle images in focus, i. e., when the light sheet thickness  $\Delta Z_0$  is smaller than the depth of field  $\delta_z$  of the optical system given by:

$$\delta_z \cong 2 f_N d_{diff} (M + 1) / M^2$$

**Depth of field:** is the region in which the image is “acceptably” sharp


M: image magnification,  $f_N$ : camera aperture,  $\lambda$ : laser wavelength

Using the diffraction limited image size of the particle

$$\delta_z \approx 4 \left(1 + 1/M\right)^2 f_N^2 \lambda$$

$$\text{For } M=0.1, f_N=5.6, \lambda=0.532 \mu\text{m} \Rightarrow \delta_z \sim 8 \text{ mm}$$

# Imaging of small particles



Since the lens is diffraction limited, the effective image diameter is larger and independent of particle size for  $d_p < 10 \mu\text{m}$


For  $d_p > 50 \mu\text{m}$  the image size is effectively the geometrical size

Actually image shape and size could differ from the Airy disk due to lens aberration

# Particles Images on CCD



- ❑ Particle images size is usually diffraction limited (Airy Disk)
- ❑ Light distribution is approximately Gaussian
- ❑ Being image sampled by a CCD sensor, Nyquist criteria requires that to be properly resolved the (diffraction limited) particle image should cover at least 4 pixel



# Particles Images on CCD


## Digitized particle image:

### Digital image

- pixel-by-pixel map of an image
- discrete 2D function





### Gray value

- value to represent brightness:  
0~255 (8-bit), 0~4095 (12-bit)



### Gray value distribution

- Particle image with Gaussian profile ( $I_p$ )
- Single pixel random noise, e.g. thermal noise ( $I_{sp}$ )
- Low frequency background noise ( $I_{lf}$ ): non-uniform illumination, flow boundary etc.
- Total intensity distribution ( $I_{tot}$ ): root-sum-square (RSS) of  $I_p$ ,  $I_{sp}$  and  $I_{lf}$



# Components of PIV System


## Light Sources

- **LASER** (stimulated light emission in an optically active media)
  - Solid state
  - Dye Laser
  - Gas
  - Semiconductor
- **GAS-DISCHARGE** (spontaneous light emission from excited gas molecule)
  - Xenon arc lamp (white light, pulsed, continuous, high power, used to simulate sunlight)
  - Deuterium arc lamp (strong UV source)
- **INCANDESCENCE** (thermal radiation)
  - Tungsten halogen lamp (visible+IR emission, high operating temperature)
- **ELECTROLUMINESCENCE** (light emission due to recombination of electron-hole pair)
  - LED** (monochromatic emission, high pulse frequency and short pulse duration, low power, no IR)
- .....

# Components of PIV System

## Light Sources – Laser Principle

- Laser is based on Light Amplification by Stimulated Emission of Radiation
- Population inversion ( $N_2 > N_1$ ) is forced to take place in a active material by means of a pump mechanism (flash lamp, electrical discharge)
- Once population inversion is achieved stimulated emission rapidly increases and an highly collimated, monochromatic light beam is emitted





# Components of PIV System

## Light Sources - Lasers

- For the illumination it is preferable to use a laser, since the laser beam is easy to form into a sheet by a cylindrical lens.
  - The laser must provide sufficient power to illuminate the seeds, and not so much that the camera becomes saturated (or damaged!).
  - This usually leads us to the use of pulsed lasers, since one obtains a high light energy during a very short time interval (typically 5 ns for a YAG-laser). Particle images will be practically frozen even for high velocities ( $> 100$  m/s).
  - If only low power is required and the camera is shuttered, a continuous laser can be used.
-

# Components of PIV System

## Light Sources - Lasers

### Continuous

Relatively low power light ( $P_{\max} = 0.01 - 50 \text{ W}$ ) of good beam quality, easy to set up and maintain

- **Helium-Neon**,  $\lambda = 633 \text{ nm}$
- **Argon Ion**, several wavelengths, 488 nm, 514 nm ..

### Single Pulse

- **Ruby**  
 $\lambda = 694.3 \text{ nm}$ , high pulse energy 1J, two pulses 1 to 500  $\mu\text{s}$  apart with one laser, 1 pulse each  $\sim 30 \text{ s}$


### Repetitive Pulse

- **Nd:YAG** (yttrium aluminum garnet,  $\text{Y}_3\text{Al}_5\text{O}_{12}$ , *heavily used in PIV, pretty cheap  $\sim 40 \text{ k€}$* )  
2<sup>nd</sup> harmonic  $\lambda = 532 \text{ nm}$ , max pulse energy  $\sim 0.5 \text{ J}$ , double pulse require 2 lasers, flash lamps 10-50 Hz, high repetition rate are also available 1-10 kHz @ 10 mJ
- **Nd:YFL** (yttrium lithium fluoride,  $\text{YLiF}_4$ ,  $\lambda = 527 \text{ nm}$ )  
high repetition rate 1-10 kHz @ 10 mJ
- **Copper Vapour**  
very expensive!, high repetition rates: 10-100 kHz, low pulse energy  $\sim 10 \text{ mJ}$

### Semiconductors, Diode-pumped Solid State Laser (DPSS)

# Components of PIV System

## Light Sources - Double Nd:YAG laser



**Fig. 2.17.** Double oscillator laser system with critical resonators

The repetition rate of a YAG-laser is typically 10-30 Hz, which is too low except for very low velocities ( $< 1$  cm/s). One therefore needs two lasers to get full freedom in terms of time separation between the pulses. Special PIV YAG-lasers are available that combine two laser cavities with a common beam outlet.

# Components of PIV System

## Light Sources – Laser Safety Issues

### ■ **SOME LASER SAFETY ISSUES**

*A typical Nd:Yag laser pulse has a energy of 0.01 to 0.5 J delivered in about 10 ns....this means 1-50 MW of laser beam peak power!!!*

- ❑ *Laser light can be VERY dangerous !!!*
- ❑ *Even laser light reflections can be VERY harmful to eyes !!!*
- ❑ *Eye damage depend upon laser power/energy and wavelength*
- ❑ *Pulsed laser light at rate below 50 Hz can be quite harmful*

### ■ **SOME VERY BASIC SUGGESTIONS**


- ❑ *Wear goggles*
- ❑ *Operate at low laser light level and bright ambiance for all alignments*
- ❑ *BE WISE!*
- ❑ *Follows regulations!*

---

*This section is not exhaustive!! it is given just to highlight safety issues!*

# Components of PIV System

## Light-Sheet Formation Optics



a) constant sheet width

b) linearly expanding sheet width

- ❑ Due to beam divergence and diameter ( $\sim 1$  cm) one single lens is not sufficient to provide a thin laser sheet when using Nd-Yag laser
- ❑ Lens configuration can be chosen by using geometrical optics laws

The Rayleigh length  $Ra$  is the length over which the light sheet has a minimum thickness  $w$

# Components of PIV System

## Image Recording Systems

- ❑ The positions of particles entrained in the flow is recorded by a camera, which is oriented 90 degrees to the plane of the light sheet.
- ❑ Early PIV was done on photographic film.
- ❑ Film still offers superior resolution and can be more sensitive but film limits the user to autocorrelation measurements, which introduces directional ambiguity.
- ❑ ***Charge coupled device (CCD) most widely used.*** Most CCD cameras frame at video rates (30 frames per second, of Hz). Often, this is not fast enough. As a result, many PIV cameras use “frame straddling”: The main challenge in frame straddling is getting the data from the first shot off the camera before the next one is taken. This is why PIV camera cost more
- ❑ Obviously pixels means resolutions. Bit depth means sensitivity and less limitation to the field of view and laser power.

---

### Important considerations in camera choice:

Spectral response, spatial resolution, temporal resolution, dynamic range, cost

# Components of PIV System

## Cross-Correlation CCD Camera

To be able to acquire two single exposed images with a time separation of the order of micro seconds, one uses a so-called *full-frame interline transfer progressive scan* CCD camera, also called a *cross-correlation* CCD-camera.

The basic idea is that the image exposed by the first laser pulse is transferred very rapidly to light-hidden areas on the CCD-chip. This is done on a pixel-by-pixel basis, i.e. each pixel has its own storage site in immediate vicinity of the light sensitive pixel area.

After the second exposure, both images are transferred to the computer.

Since a lot of data has to be transferred, it is only possible to take a few double images per second.

The temporal resolution of the flow is thus in general very poor with this technique.


---



# Components of PIV System

## Cross-Correlation CCD Camera

- **CCD chip contains light sensitive cells and storage cells in an  $n$  by  $m$  array**
- **After the first exposure the image information is transferred to the storage cell and the CCD chip is ready for the second exposure**
- **Minimum time between cross correlation exposures is 200 ns.**



Electrons generated by impinging photon are continuously stored in the pixel site until read out

# Components of PIV System

## PIV with Cross-Correlation CCD Camera

**First laser pulse**

Duration ~5 ns

**Second laser pulse**

Duration ~5 ns

on

off

*Laser Pulse*

$2 / \text{frame rate}$

on

off

*Electronic Shutter*

**Photodiode array  
is charged**

**Overexposure  
due to register  
delay**  
*About 30-80 ms*


**Register clears;  
Second image is transferred**

**Charges moved from  
pixel to adjacent  
storage cell**

Duration ~ 0.3  $\mu\text{s}$







**Second Image is captured,**  
but electronic shutter must stay open until  
register clear

**First image is transferred**



# Components of PIV System

## PIV with CCD Camera (general)



PIV frame rate is (camera frame rate)/2

# Components of PIV System

## Camera Optics


Typically, SLR lenses are used, and these also impact the results. Good quality lenses are mandatory to obtain accurate measurements.

The field of view and the working distance are functions of the camera optics.

Usually a narrow band filter is used in front of the lenses (or between lenses and camera) to block ambient light.



# The Principle of PIV



# The Principle of PIV




## Single Exposure – Double Frame

Two light pulses separated by short time interval

At each pulse the positions of particles are recorded using a digital or film camera (one pulse in one frame).

Two images are needed to evaluate the particle velocity.

Too long illumination times will result into streaks of the particle images that will not allow to determine the exact particle location in the fluid.



# The Principle of PIV




## Multiple Exposure – Single Frame

Two or more light pulses separated by a short time interval

At each light pulse the positions of particles are recorded using a digital or film camera (all pulses on a single frame).


Too long illumination times will result into streaks of the particle images that will not allow to determine the exact particle location in the fluid.





# Analysis of PIV images

## Optical set-up



Particle image velocimetry (PIV) allows the mean velocity of a fluid to be measured in a small area  $D_I$ , illuminated by a light sheet and seeded with tracer particles of concentration  $C$




# Analysis of PIV images

PIV images are defined accordingly to their particle density


Particle density is represented by two dimensionless parameters:

- Image particle density  $N_I$  represents the number of individual particles in the “probe volume”.
- Source density  $N_S$  represents whether the particle images are overlapping ( $N_S > 1$ ) or can be recognized individually ( $N_S < 1$ ).

$$N_S = \frac{C \Delta z_0}{M_0^2} \cdot \frac{\pi}{4} d_\tau^2$$




$$N_I = \frac{C \Delta z_0}{M_0^2} \cdot D_I^2$$



# Analysis of PIV images

Particle density defines the type of velocimetry that can be used:

- ❑ Particle Tracking Velocimetry (PTV)
- ❑ Particle Image Velocimetry (PIV).



Respect to low-image density PIV the interrogation in any location normally contains sufficient particle images to obtain a valid measurement.


In general, the number of velocity data extracted from high-image-density PIV is an order of magnitude larger than for low-image-density PIV

# Analysis of PIV images




Several recording strategies can be implemented


single frame  
double exposure



single frame  
multiple exposure



multiple frame  
single exposure




Since it is impossible to track the individual particles, the *mean* displacement of particles in a small region of the image(s) (the interrogation area) has to be calculated by a numeric algorithm, es. spatial correlation in case of digital PIV (DPIV).

---


# Analysis of PIV images

## Double exp-single frame (Autocorrelation)



**Double-exposed  
image**

**Interrogation  
region**



**Spatial  
correlation**

# Analysis of PIV images

## Double exp-single frame (Autocorrelation)

For a double exposed image, the auto-correlation algorithm is applied to a small interrogation area ( $D_I$ ) of size  $M \times N$ . the **auto-correlation** function is given as


$$\Phi(m, n) = \sum_{i=1}^M \sum_{j=1}^N I_1(i, j) \cdot I_1(i + m, j + n)$$

where  $I_1$  is the grey value distributions

Auto correlation show **one peak at zero displacement** and two mirrored displacement peaks.

The mean particle image displacement in the interrogation area is determined by the position of the max value of  $\Phi(m, n)$ .

Auto correlation cannot resolve directional ambiguity by itself.



Autocorrelation field

# Analysis of PIV images


## Single exp–double frame (Cross-correlation)

Cross-correlation algorithm using two images of particles in the flow is used for calculating the velocity vector

The 2-D average particle motion is obtained by correlating the interrogation area's of the two images.

$$\Phi(m, n) = \sum_{i=1}^M \sum_{j=1}^N I_1(i, j) \cdot I_2(i + m, j + n)$$


-  $m$  and  $n$  are shifts of the second IA respect to the first one





# Analysis of PIV images

## Single exp–double frame (Cross-correlation)




Average particles displacement is given by position of the peak

The cross-correlation method contains two assumptions:

- 1) The contribution from noise is negligible
- 2) The particles within the interrogation area move uniformly in a single direction during their double image capturing time


Since velocity gradients may prevent the validity of the 2nd assumption, a condition on the velocity variation is recommended



CROSS-CORRELATION PROCEDURE

# Analysis of PIV images

## Single exp–double frame (Cross-correlation)



For a single exposed PIV recording pair, the cross-correlation algorithm is applied to a small interrogation area ( $D_1$ ) of size  $M \times N$ . Two samples are chosen, from the first and second images and the *cross-correlation* function is given as

$$\Phi(m, n) = \sum_{i=1}^M \sum_{j=1}^N I_1(i, j) \cdot I_2(i + m, j + n)$$

where  $I_1$  and  $I_2$  are grey value distributions

The mean particle image displacement in the interrogation area is determined by the position of the max value of  $\Phi$ .

Cross correlation renders only the velocity peak.

Cross correlation resolves directional ambiguity by using two consecutive images.




# Analysis of PIV images

## Single exp–double frame (Cross-correlation)

Mean Intensity    Fluctuating    Displacement


↓                      ↓                      ↙

$$\Phi = R(s) = R_C(s) + R_F(s) + R_D(s)$$



# Analysis of PIV images

## Single exp–double frame (Cross-correlation)



**Fig. 5.102a–c** Separation of the displacement-correlation term  $R_D$  (a) into mean (b) and fluctuating terms (c)

$$R_D(s) = \underbrace{\langle R_D(s) | u \rangle}_{\text{Displacement}} + \underbrace{R'_D(s)}_{\text{Noise}}$$

The correct displacement-correlation peak can be identified only when its amplitude is larger than the highest random correlation peak in  $R'_D(s)$ .

Otherwise the interrogation analysis yields a *spurious* value for the particle-image displacement.

The ratio of the highest correlation peak and that of the second highest correlation peak is a measure of the *detectability* of the displacement-correlation peak, indicated by  $D_0$ , and can be considered as a lower limit of the signal-to-noise ratio.

# Analysis of PIV images


## Particle Pair Loss

### IN-PLANE LOSS ( $F_I$ )



Particle pair are lost due to particle in-plane displacement to laser sheet.

### OUT-PLANE LOSS ( $F_O$ )


Particle pair are lost due to particle displacement perpendicular to laser sheet.



**Fig. 5.104** The spatial correlation for increasing in-plane displacement (at constant image density)




**Fig. 5.105** The spatial correlation for increasing out-of-plane motion



# Analysis of PIV images

## Correlation based techniques




$N_1$  = image density  
 $F_1$  = in plane loss of images  
 $F_0$  = out of plane loss of images

If the value of  $N_1 F_1 F_0$  is above 5 the probability to make a valid measurement is above ~90%

**Fig. 5.107** The valid detection probability for the displacement-correlation peak as a function of the image density  $N_1$  and the in-plane and out-of-plane loss of correlation,  $F_0$  and  $F_1$ , respectively (after Keane and Adrian [5.434])

# Analysis of PIV images

## Sub-Pixel Resolution




- ❑ The displacement obtained by cross-correlation (auto-correlation) algorithms is resolved within  $\pm 0.5$  pixel (error could be as large as 1 pixel).
- ❑ Under this circumstance the relative error for digital PIV measurements would not be better than  $1/8$  ( $\approx 13\%$ ) when using an IA of  $32 \times 32$  pixels and a displacement of 8 pixels
- ❑ A sub-pixel resolution can be achieved by interpolating the discrete cross-correlation (auto-correlation) function (uncertainty  $\sim 0.05$ – $0.10$  pixel)
- ❑ Intensity distribution of particle image has a Gaussian shape, a Gaussian fit is applied to the 3 adjacent points of  $\Phi$  (along  $x$  and  $y$ )




# Analysis of PIV images

## Sub-Pixel Resolution



Gaussian Intensity distribution



**Fig. 5.108** Detailed shape of the correlation peak (in one dimension) as a function of the subpixel part of the displacement (here shown for a 1.6 pixel particle-image diameter). The *horizontal shaded region* represents the background noise level in the correlation

$$\hat{\varepsilon}_X = \frac{\ln R_{-1} - \ln R_{+1}}{2(\ln R_{-1} + \ln R_{+1} - 2 \ln R_0)} ,$$

with

$$R_{-1} \equiv R_D[m_0 - 1, n_0], \quad R_0 \equiv R_D[m_0, n_0],$$

$$R_{+1} \equiv R_D[m_0 + 1, n_0] ,$$

$$I(X, Y; t) = \frac{2\alpha}{\frac{\pi}{4} d_\tau^2} \exp \left[ -8 \frac{(X - X_0)^2 + (Y - Y_0)^2}{d_\tau^2} \right]$$

- ❑ A sub-pixel resolution can be achieved by interpolating the discrete cross-correlation (auto-correlation) function (uncertainty ~0.05–0.10 pixel)
- ❑ Intensity distribution of particle image has a Gaussian shape, a Gaussian fit is applied to the 3 adjacent points of F (along x and y)

# Analysis of PIV images

## Correlation based techniques FFT


- ❑ The computation of correlation is very expensive

number of operation  $\sim (M \times N)^2 \sim 10^6$

- ❑ Fast Fourier Transform (FFT) algorithm allows to speed up the correlation process ( $4MN \log_2 MN$ ,  $\sim 10^4$ )


- ❑ FFT could introduce artifact in  $\Phi(m, n)$

$$\Phi(m, n) = \sum_{i=1}^M \sum_{j=1}^N I_1(i, j) \cdot I_2(i + m, j + n)$$



# Analysis of PIV images

## Single exp–double frame (Cross-correlation)





# Performance of Correlation Techniques

## Random Errors




Fig. 1. The RMS estimation error for the displacement as a function of the displacement  $u$  in pixel units for a  $32 \times 32$ -pixel interrogation region in digital PIV cross-correlation analysis, for particle-images with a diameter of 2 and 4 pixels with an image density of 10, with zero out-of-plane displacement. The solid lines are analytical results (Westerweel 1993a, b); the symbols are obtained from Monte Carlo simulations (the error bar represents the uncertainty of the simulation result)

Random Errors are approximately independent of particle displacement


At displacement below  $<0.5$  pixels random error decrease linearly

Random error is proportional to particle image diameter

$$\sigma_{\Delta X} \cong c d_\tau$$


# Performance of Correlation Techniques

## PIV Peak Locking (Bias Error)




For particle size less than 2 pixel the displacement is biased towards integer values

The three point gaussian peak-fit is unsuited at these particle image diameter



# Performance of Correlation Techniques

## PIV Bias Errors




Estimated displacement are biased towards smaller value

- ❑ In plane loss of pare (large displacements gives less strong signal)
- ❑ Finite size of interrogation area

Divide correlation with factor found from correlation of uniform signal


# Performance of Correlation Techniques

## Optimization of Particle Image Diameter



**Fig. 5.112** The rms displacement error as a function of the particle image diameter (in pixel units) for  $32 \times 32$  px interrogation windows and constant image density  $N_I = 15$  (open symbols) and constant source density  $N_s = 0.05$ . The solid line represents a theoretical result, while the dotted and dashed lines represent scaling relations. (After Westerweel [5.433])

For cross-correlation technique and three point Gaussian peak-fit the optimum particle size is about 2-4 pixels



# Performance of Correlation Techniques

## Interrogation Area Optimization

To minimize loss of pair and bias errors the size of Interrogation Area should fulfill the following conditions

### *Interrogation Area optimization*

- Average number of particle images

(C: number of particles per unit volume)

$$N = C \Delta z \pi [D_I / (2M)]^2 > 5 \div 10$$

- Interrogation Area dimension:

$$(u^2 + v^2)^{1/2} \Delta t < D_I / 4M$$

- Out-of-plane velocity component

( $\Delta z$ : laser sheet thickness)

$$w \Delta t < \Delta z / 4$$

- Spatial gradients

$$|\Delta u| / |u| < 0.2$$

In high *image density* ( $N_I$ ) there are less problems than in low image density cases, where the experimental parameters need to be optimized.

# Performance of Correlation Techniques

## Dynamic Velocity/Spatial Range

dynamic spatial range  $DSR = \frac{\text{image format}}{\text{interrogation resolution}}$   $\rightarrow L_x / \Delta X_{\max} \sim L_x / D_I$

dynamic velocity range  $DVR = \frac{\text{max. particle displacement}}{\text{displacement error}}$   $\rightarrow \frac{\Delta X_{\max} - \Delta X_{\min}}{\Delta X_{\min}} = \frac{D_I}{4\Delta X_{\min}} - 1.$

$\Delta X_{\min} = 0.1 \text{ px}$

system	format	DSR	DVR
standard video	768×494	50	140
large-format video	2029×1044	130	140
35 mm film	7000×4800	440	140
4"×5" film (Tmax)	25000×20000	1560	140
8"×10" film (TP)	160000×128000	10000	140

# Advanced Digital Interrogation Technique




## Velocity maps of a swirling flow

Recursive correlation processing allows to iteratively improve spatial resolution.

Recursive correlation processing uses first step correlation to obtain a mean displacement and direction. These data are used to define the centres of a search window for image pairing.

By offsetting the second interrogation windows the matched particle image increase (SNR of correlation peak increase).

Thus size of interrogation area can be reduced at each step without significantly affect the measurements uncertainty






# Advanced Digital Interrogation Technique



The pictures illustrate the possibility of processing PIV images to the limits of optical resolution, even when significant out of plane flow exists






# Developments of PIV

## 3-D PIV

Stereoscopic planar PIV measures the 3 dimensional velocity field in a plane


Sets of vector maps are recorded and calculated by cross correlation

The 3D vectors are reconstructed from projections of the two maps to the light sheet plane



# Developments of PIV

## 3-D PIV - Fundamentals of stereo vision



True 3D displacement is estimated from a pair of 2D displacements as seen from left and right camera respectively


# Developments of PIV

## 3-D PIV

Focusing an off-axis camera requires tilting of the CCD-chip (Scheimpflug condition)


3D evaluation requires a numerical model, describing how objects in space are mapped onto the CCD-chip of each camera

Parameters for the numerical model are determined through camera calibration



# Developments of PIV

## 3-D PIV – Calibration



Images of a calibration target are recorded.

The target contains calibration markers in known positions.

Comparing known marker positions with corresponding marker positions on each camera image, model parameters are adjusted to give the best possible fit.

---


# Developments of PIV

## 3-D PIV – Overlapping fields of view

3D evaluation is possible only within the area covered by both cameras.

Due to perspective distortion each camera covers a trapezoidal region of the light sheet.

Careful alignment is required to maximize the overlap area.



# Developments of PIV

## Time-Resolved PIV

Same principle as classical PIV


High frame rate with special laser and camera

Time resolved flow analysis

### Hardware

Cw-diode pumped Nd:Yag laser

High-speed megapixel CMOS sensor



# Developments of PIV

## Micro-PIV

High Spatial resolution 100-10  $\mu\text{m}$

Particle smaller than laser light wavelength (Fluorescent seeding)

Volume illumination

Low particle density

Special correlation technique




Fig. 1 Diagram of typical micro-PIV system.

# PIV Applications


PIV technique has been applied to investigate a very broad range of technical and scientific problems.

- ❑ Gasdynamics, Fluidynamics, Combustion & Turbulence
  - Flow across shock waves
  - Turbomachinery (flow around turbine and compressor blades)
  - Laminar and turbulent flames
  - Swirl combustor (aerospace propulsion and stationary power systems)
  - Premixed and non premixed burners
- ❑ Aerospace
  - Flow inside rocket motors
  - Flow over wings
  - Wake of hovering rotor
- ❑ Car industry
  - Vehicle Aerodynamics
- ❑ Naval Hydrodynamics
  - Ship wakes
  - Propeller flows
  - Underwater vehicle hydrodynamics
- ❑ Bioengineering
  - Flow in artificial heart valves
  - Flow in lung
- ❑ Appliance
- ❑ .....




# Examples of PIV Measurements

## Measurements In a Swirling Air Flow



### Burner geometry




### Burner - top view



- Swirl generated by axial + tangential air injection
- Transverse fuel mixture injection
- Fuel injection below the burner exit (partial air/fuel mixing)

# Examples of PIV Measurements

## Measurements In a Swirling Air Flow



2-D mean velocity pattern in the axial section of a axis-symmetric swirling air flow.

Burner with two co-axial flows: inner jet and co-flow swirling air

A large toroidal recirculation region is generated in the central region.


Average over 200 double images

Green = positive velocities

Bleu= negative velocities

# Examples of PIV Measurements

## Measurements In a Swirling Air Flow




2D plot of the axial velocity component measured by PIV in a swirled natural gas flame (50%  $H_2$ ) and streamlines



# Examples of PIV Measurements



## Average and instantaneous vortex core in a Swirling Air Flow




PIV measurements made in a swirling jet flow under PVC regime

# Examples of PIV Measurements

## Transonic Wind Tunnel ONERA - Modane



*Figure 1 : Profil d'aile dans la soufflerie transsonique S2 Modane.*



*Figure 2 : Dispositif de génération de plan laser.*

# Examples of PIV Measurements



## Transonic Wind Tunnel - ONERA - Modane





Figure 3 : Cartes de vitesses moyennes à Mach 0,82 et incidence 3.7° et 4.5°.

Mean velocity distribution in a longitudinal plane of the wing extrados. PIV measurements show a shock on the wing extrados

# Examples of PIV Measurements

## 3-D PIV in the wake of a hovering rotor





Figure 3: Rotor test stand and PIV apparatus.

# Examples of PIV Measurements

## 3-D PIV in the wake of a hovering rotor





# References



- ❑ Raffel, M. et al. *Particle Image Velocimetry. A Practical Guide*, Springer 2007
- ❑ *Handbook of Experimental Fluid Mechanics*, edited by Tropea C., Yarin A.L., Foss J.F., Springer 2007
- ❑ *Particle Image Velocimetry, New Developments and Recent Applications* Series: Topics in Applied Physics , Vol. 112, Schröder, Andreas; Willert, Christian E. (Eds.) 2008
- ❑ Adrian R.J. *Twenty years of particle image velocimetry*, Experiments in Fluids, 2005
- ❑ Grant I. *Particle image velocimetry: a review*, IMechE 1997